

Stools and Chairs



Distributore Ufficiale per l'Italia

REF ITALIA S.R.L

Via Nomentana, 403 – 00162 ROMA

Tel. +39 0645547880 – Fax +39 0645547881

Mail: info@refitalia.it

www.refitalia.it – www.colposcopi.it





Stools and Chairs

Swivel Stools

The swivel stools belong to the indispensable equipment in every medical practice or facility. To ensure users and patients enjoy a comfortable sitting position, these swivel stools can be continuously adjusted in height (by foot or hand release).



1 DH-2000 (soft double castors)



2 DH-2000/RA (soft double castors)



3 DHMF-2000/LG



4 DHFA-2000

Article Number	DH-2000	DH-2000/RA	DH-2000/LG	DH-2000/RALG	DHMF-2000/LG	DHMF-2000/RALG	DHFA-2000
Picture(s)	1	2	5	Similar 2	3	Similar 6	4
Price							
Maximum load (kg)	140		140		140		140
Upholstery thickness - Seat (mm)	75		75		75		75
Seat diameter (mm)	380		380		380		380
Seat height adjustment (mm)	430-550		520-700		520-700		460-580
Height adjustment by							
▪ Foot release	–	–	–	–	–	–	✓
▪ Hand lever release	✓	–	✓	–	✓	–	–
▪ Ring release (under the seat)	–	✓	–	✓	–	✓	–
Height-adjustable foot ring (Ø 455 mm)	–		–		✓		–
Aluminium base - Polished (five legs, Ø 610 mm)	✓		✓		✓		–
Aluminium base - Black (five legs, Ø 610 mm)	DH-1000	DH-1000/RA	DH-1000/LG	DH-1000/RALG	DHMF-1000/LG	DHMF-1000/RALG	✓

✓ included in delivery ▪ available (extra charge) – not available

Additional Features

- The swivel stools are equipped as a standard with braked (hard) double castors (Ø 50 mm) for soft floors. They are also available with soft castors (for hard floors) at no extra charge. Please specify upon ordering.
- The swivel stools may be ordered with gliders (feet), at no extra charge. Please note that, with gliders, the total height decreases by approximately 25 mm.
- On the swivel stool DHFA-2000 with foot ring release, the height adjustment of the seat cushion is activated with the foot via a chrome-plated foot ring (Ø 300 mm), located under the base, between the feet of the stool.
- On all other swivel stools, the height adjustment is done by hand (with a ring release or a hand lever).

Options

- The OR swivel stools (with electrically conductive double castors and electrically conductive synthetic leather upholstery) are shown on page 20.



5 DH-2000/LG



6 DHMF-1000/RALG

Stools and Chairs

Swivel Chairs

The swivel chairs are a further development of the swivel stools presented on the previous page. They are equipped with a height and tilt-adjustable backrest. To ensure users or patients enjoy a comfortable sitting position, the swivel stools can also be continuously adjusted in height (by foot or hand release).



1 DS-2000/RL



2 DS-2000/RLLG (soft double castors)



3 DSMF-2000/RLLG



4 DSFA-2000/RL

Article Number	DS-2000/RL	DS-2000/RLLG	DSMF-2000/RLLG	DSFA-2000/RL
Picture(s)	1	2	3	4
Price				
Maximum load (kg)	140	140	140	140
Upholstery thickness - Seat (mm)	75	75	75	75
Upholstery thickness - Back (mm)	45	45	45	45
Seat diameter (mm)	380	380	380	380
Seat height adjustment (mm)	460-580	550-730	550-730	490-610
Height adjustment mechanism				
▪ Foot release	–	–	✓	✓
▪ Hand lever release	✓	✓	–	–
Height-adjustable foot ring (Ø 455 mm)	–	–	✓	–
Aluminium base - Polished (five legs, Ø 610 mm)	✓	✓	✓	✓
Aluminium base - Black (five legs, Ø 610 mm)	DS-1000/RL	DS-1000/RLLG	DSMF-1000/RLLG	–

✓ included in delivery ▪ available (extra charge) – not available

Additional Features

- These swivel chairs are equipped as a standard with braked (hard) double castors (Ø 50 mm) for soft floors. They are also available with soft castors (for hard floors) at no extra charge. Please specify upon ordering.
- Additionally, they may be ordered with gliders (feet), at no extra charge. Please note that, with gliders, the

total height decreases by approximately 25 mm.

- On model DSFA-2000/RL with foot ring release, the height adjustment is done with the foot via a chrome-plated foot ring (Ø 300 mm), located under the base, between the feet of the chair. On all other swivel chairs, the height adjustment is done with a hand lever.

- The seat has a slight inclination of 3° towards the backrest.

Options

- The OR swivel chairs (with electrically conductive castors and synthetic leather upholstery) are shown on catalogue page 20.



5 DS-1000/RLLG

Stools and Chairs

Office Chairs

The office chairs are similar to our blood collection chairs (catalogue page 101). They also feature an ergonomic backrest which is continuously adjustable in height and tilt and an ergonomic seat for added comfort. Furthermore, they are available with or without armrests.



1 BR-2000



2 BR-1000



3 B-2000



4 B-1000

Article Number	BR-2000	BR-1000	B-2000	B-1000
Picture(s)	1	2	3	4
Price				
Maximum load (kg)	140		140	
Aluminium base (five legs, Ø 610 mm)	Polished	Black	Polished	Black
Upholstery thickness - Seat (mm)	75		75	
Seat dimensions (mm)	460x450 (WxD)		460x450 (WxD)	
Upholstery thickness - Back (mm)	40		40	
Dimensions - Back (mm)	450x500 (WxH)		450x500 (WxH)	
Height adjustment by hand lever (mm)	420-540 (Gas pressure spring)		420-540 (Gas pressure spring)	
Pair armrests	✓		—	
Set of braked hard double castors (Ø 50 mm)	✓		✓	
Gliders (height 420-540 mm)	SGL (No extra charge)			
Height adjustment 510-690 mm	LG			
Height-adjustable foot ring	MF			

✓ included in delivery ▪ available (extra charge) — not available

Options

- The office chairs are equipped as a standard with braked (hard) double castors (Ø 50 mm) for soft floors.
- They are also available with soft castors (for hard floors) at no extra charge. Please specify upon ordering.



5 Base with gliders

Stools and Chairs

Examination Stools with Saddle Seat and Square Stools

The examination stools with saddle seat are based on our swivel stools (page 104). They differ only in the black, PU-foamed, ergonomically shaped, adjustable saddle seat. The advantage over a standard seat, due to the individual adjustment of the pelvic tilt and the associated automatic straightening of the upper body (spine), is that the posture can be optimized and back tension thus be minimized.

The square examination stools are an alternative to the various stools and chairs described previously. They are available in a lacquered or chrome-plated frame finish.



1 SH-2000/LG (soft double castors)



2 SH-1000/LG



3 v-2000



4 v-1000

Article Number	SH-2000	SH-1000	SH-2000/LG	SH-1000/LG	V-2000	V-1000
Picture(s)	Similar 1	Similar 2	1	2	3	4
Price						
Maximum load (kg)	140		140		130	
Square steel tubing profiles, 25x25 / 20x20 mm	-		-		Chrome	Lacquered
Aluminium base (five legs, Ø 610 mm)	Polished	Black	Polished	Black	-	
Upholstery thickness - Seat (mm)	-		-		67	
Seat dimensions (mm)	420x370 (WxD)		420x370 (WxD)		400x400 (WxD)	
Seat height / Height Adjustment by hand lever (mm)	460-580 (Gas pressure spring)		550-730 (Gas pressure spring)		500	
Set of braked hard double castors (Ø 50 mm)	✓		✓		-	
Gliders (feet)	SGL (No extra charge)		SGL (No extra charge)		-	
Height adjustment range (mm)	435-555		525-705		-	

✓ included in delivery ▪ available (extra charge) - not available

Options

- The saddle seat stools are equipped as a standard with braked (hard) double castors (Ø 50 mm) for soft floors.
- They are also available with soft castors (for hard floors) at no extra charge, please specify upon ordering.