

Sleek compact design fits into any environment



Tilt-and-swivel monitor for optimal viewing



Customizable control panel to personalize your user interface



Completely sealed control panel to simplify clean-up



A Clear Vision

Ultrasound reimaged for ease of use

Born of a vision to deliver meaningful design innovations that benefit the user, the Acclarix AX4 features a host of design breakthroughs that make day-to-day operation in general imaging environments easy, fast and intuitive. The result is an elegant simplicity where form and function meet at the tips of your fingers.

- Best-in-class image quality with one-key optimization
- Application-specific exam and imaging presets
- Unique gesture-control user interface
- Tilt-and-swivel HD monitor
- Completely sealed control panel aids in maintaining infection control
- Needle visualization technology improves needle identification

Transducers and Options

<p>L12-5Q</p>  <p>5 - 12 MHz Linear Array</p>	<p>L17-7SQ</p>  <p>7 - 17 MHz High Frequency Compact Linear Array</p>	<p>C5-2Q</p>  <p>2 - 5 MHz Curved Linear Array</p>	<p>C5-2MQ</p>  <p>2 - 5 MHz 3D/4D Mechanical Curved Array</p>
<p>E8-4Q</p>  <p>4 - 8 MHz Endocavity Curved Array</p>	<p>MCS-4Q</p>  <p>4 - 8 MHz Micro Convex Array</p>	<p>P5-1XQ</p>  <p>1 - 5 MHz Phased Array</p>	<p>P5-1Q</p>  <p>1-5 MHz Phased Array</p>

 <p>Height-adjustable 4-wheel swivel cart</p>	 <p>Multi-transducer connector option</p>	 <p>Sturdy, versatile travel case</p>
---	---	---

Redefining Value for General Imaging Ultrasound

Acclarix AX4
Compact Ultrasound System




U.S. inquiries:
Edan Diagnostics, Inc. | 9833 Pacific Heights Blvd., Suite E/F
San Diego, CA 92121-4707 | +1.858.750.3066
www.edandiagnostics.com | edan-info@edandiagnostics.com

Inquiries outside of the U.S.:
Edan Instruments, Inc. | No.15 Jinhui Rd., Jinsha Community,
Kengzi Subdistrict, Pingshan District, Shenzhen | 518067 P.R. China
+86.755.26898326 | www.edan.com.cn | info@edan.com.cn




© 2017 Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice.







This fluid collection on the anterolateral aspect of the knee is clearly defined with the L12-5Q transducer.



TAI for color Doppler enhances flow sensitivity and spatial resolution as demonstrated with this thyroid.



The C5-2Q transducer delivers exceptional color and spectral Doppler sensitivity visualizing the left hepatic vein.



The 3D/4D feature enables fast and easy acquisition for outstanding fetal face detail.

Stunning Clarity

Definitive imaging

Sleek design and user-friendly features are just part of the story. The Acclarix AX4 is designed to deliver amazing clarity and stunning image quality quickly—eliminating the need to tweak imaging controls.

- High fidelity architecture results in superb resolution
- Tissue Adaptive Imaging (TAI) continuously and automatically optimizes imaging allowing more focus on the patient
 - In B-mode, TAI fine tunes multiple parameters to provide the best possible image quality
 - In Doppler, TAI automatically adjusts for flow state providing improved continuity, border detection and fill-in

Redefining innovation through value and performance

Introducing a new price performance standard for ultrasound. The remarkable Acclarix™ AX4 Compact Ultrasound System delivers a powerhouse combination of imaging performance and ease-of-use to meet the specific demands of general imaging ultrasound. The Acclarix AX4 has been designed from the ground up with a relentless focus on delivering unexpected levels of performance at a very attractive price point.

- Definitive Imaging
- Distinctive Design
- Intelligent Workflow
- Intrinsic Quality



The Virtue of Value

How can Edan deliver so much innovation AND so much value? By capitalizing on the experience and expertise of a truly global enterprise. Acclarix technology is the result of Silicon Valley innovation and engineering excellence combined with Chinese design and manufacturing proficiency. The Acclarix AX4 delivers unmatched value and performance across a broad range of applications.

- Abdomen
- Breast
- Prostate
- Difficult-to Image
- OB/GYN
- Nerve
- Vascular
- Musculoskeletal
- Cardiac Screening
- Small Parts
- Intraoperative